

Bearing Sensor bracket

Copyright

Ergonomic positioning!



New mobile support has been specially designed to place a digital sensor in any position.

Ideal for emergency services, you can place it under a table or a stretcher, in profile, and above a bed for the pictures in your room.

Three levels of height adjustments, assisted by a hydraulic jack, you can set it to any height between 61 cm and 158 cm.

Orientation: 0, 30, 45 and 90 °.

Very compact and takes up very little space in the room, 0.33 M².

The jaws are adapted according to the dimensions of voting sensor.



Bed X-ray



Front X-ray



Profile X-ray



Chest X-ray



Members X-ray

